

ICP-MS Internal Standard for Method 200.8, 100 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-ICP-MS-200.8-IS-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Scandium (Sc)	[7440-20-2]	100	µg/mL		
	Yttrium (Y)	[7440-65-5]	100	µg/mL		
	Indium (In)	[7440-74-6]	100	µg/mL		
	Terbium (Tb)	[7440-27-9]	100	µg/mL		
	Bismuth (Bi)	[7440-69-9]	100	µg/mL		