

**ICP-MS Tuning solution for Method 200.8, 10 µg/mL in 2-5% Nitric Acid**

Art. ID ACCU-ICP-MS-200.8-TUN-1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	10	µg/mL		
	Magnesium (Mg)	[7439-95-4]	10	µg/mL		
	Cobalt (Co)	[7440-48-4]	10	µg/mL		
	Indium (In)	[7440-74-6]	10	µg/mL		
	Lead (Pb)	[7439-92-1]	10	µg/mL		