

Density Standard 0.6960g/ml @15°C

| | |
|-----------------|--|
| Art. ID | REA-DEN15010 |
| Unit | 100 mL |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2 |

Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Traceable to a National or International standard (eg OIML)

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 15°C | Density | | 0,696 | g/mL | | |