

Lead ICP-MS Standard, 100 µg/mL in diluted Nitric Acid

Art. ID ACCU-ICP-MS-29N-0.01X-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Lead (Pb) | [7439-92-1] | 100 | µg/mL | | |