

Silver ICP-MS Standard, 1000 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-ICP-MS-53N-0.1X-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 1000 | µg/mL | | |