

Density Standard 0.8428g/ml @15°C (Mixture of toluene and bihexyl)

Art. ID REA-DEN15100
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1294 II E2

Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 15°C	Density		0,8428	g/mL		