

**Density Standard 0.8428g/ml @15°C (Mixture of toluene and bihexyl)**

Art. ID REA-DEN15100  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1294 II E2

**Description**

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 15°C	Density		0,8428	g/mL		