

## **ICP-MS Interference Check Standard 1 for Method 6020, Varied conc, in 2-5% Nitric Acid**

Art. ID ACCU-ICP-MS-6020-INT1-1  
Unit 100 ML  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	1000	µg/mL		
	Chloride (Cl-)		10000	µg/mL		
	Calcium (Ca)	[7440-70-2]	1000	µg/mL		
	Carbon (C)	[7440-44-0]	2000	µg/mL		
	Iron (Fe)	[7439-89-6]	1000	µg/mL		
	Magnesium (Mg)	[7439-95-4]	1000	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	20	µg/mL		
	Phosphorus (P)	[7723-14-0]	1000	µg/mL		
	Potassium (K)	[7440-09-7]	1000	µg/mL		
	Sodium (Na)	[7440-23-5]	1000	µg/mL		
	Sulfur	[7704-34-9]	1000	µg/mL		
	Titanium (Ti)	[7440-32-6]	20	µg/mL		