

ICP-MS Tuning Solution for Method 6020, 10 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-ICP-MS-6020-TUN-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cobalt (Co)	[7440-48-4]	10	µg/mL		
	Indium (In)	[7440-74-6]	10	µg/mL		
	Lithium (Li)	[7439-93-2]	10	µg/mL		
	Thallium (Tl)	[7440-28-0]	10	µg/mL		