

Tungsten ICP-MS Standard, 1000 µg/mL in Water, tr Ammonium hydroxide

Art. ID ACCU-ICP-MS-65W-0.1X-1
Unit 100 mL
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Tungsten (W) | [7440-33-7] | 1000 | µg/mL | | |