

Density Standard 1.0337g/ml @15°C

Art. ID	REA-DEN15150
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1105 III E1

Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 15°C	Density		1,0337	g/mL		