

ICP-MS Quality Control Sample 2, 10 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-ICP-MS-QC2-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10	µg/mL		
	Antimony (Sb)	[7440-36-0]	10	µg/mL		
	Arsenic (As)	[7440-38-2]	10	µg/mL		
	Barium (Ba)	[7440-39-3]	10	µg/mL		
	Beryllium (Be)	[7440-41-7]	10	µg/mL		
	Cadmium (Cd)	[7440-43-9]	10	µg/mL		
	Calcium (Ca)	[7440-70-2]	10	µg/mL		
	Chromium (Cr)	[7440-47-3]	10	µg/mL		
	Cobalt (Co)	[7440-48-4]	10	µg/mL		
	Copper (Cu)	[7440-50-8]	10	µg/mL		
	Iron (Fe)	[7439-89-6]	10	µg/mL		
	Lead (Pb)	[7439-92-1]	10	µg/mL		
	Magnesium (Mg)	[7439-95-4]	10	µg/mL		
	Manganese (Mn)	[7439-96-5]	10	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	10	µg/mL		
	Nickel (Ni)	[7440-02-0]	10	µg/mL		
	Potassium (K)	[7440-09-7]	10	µg/mL		
	Selenium (Se)	[7782-49-2]	10	µg/mL		
	Silver (Ag)	[7440-22-4]	10	µg/mL		
	Sodium (Na)	[7440-23-5]	10	µg/mL		
	Thallium (Tl)	[7440-28-0]	10	µg/mL		
	Thorium (Th)	[7440-29-1]	10	µg/mL		
	Uranium (U)	[7440-61-1]	10	µg/mL		
	Vanadium (V)	[7440-62-2]	10	µg/mL		
	Zinc (Zn)	[7440-66-6]	10	µg/mL		