

ICP-MS Spiking Solution Set

Art. ID	ACCU-ICP-MS-SPIKE-1-SET
Unit	2 x 100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG)

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
ACCU-ICP-MS-SPIKE-S-1	Antimony (Sb)	[7440-36-0]	100	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Arsenic (As)	[7440-38-2]	50	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Barium (Ba)	[7440-39-3]	250	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Beryllium (Be)	[7440-41-7]	25	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Cadmium (Cd)	[7440-43-9]	50	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Chromium (Cr)	[7440-47-3]	250	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Cobalt (Co)	[7440-48-4]	100	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Copper (Cu)	[7440-50-8]	250	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Lead (Pb)	[7439-92-1]	100	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Nickel (Ni)	[7440-02-0]	125	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Selenium (Se)	[7782-49-2]	25	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Silver (Ag)	[7440-22-4]	25	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Thallium (Tl)	[7440-28-0]	25	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Vanadium (V)	[7440-62-2]	150	µg/mL		
ACCU-ICP-MS-SPIKE-S-1	Zinc (Zn)	[7440-66-6]	250	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Antimony (Sb)	[7440-36-0]	100	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Arsenic (As)	[7440-38-2]	50	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Barium (Ba)	[7440-39-3]	250	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Beryllium (Be)	[7440-41-7]	25	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Cadmium (Cd)	[7440-43-9]	25	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Chromium (Cr)	[7440-47-3]	100	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Cobalt (Co)	[7440-48-4]	100	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Copper (Cu)	[7440-50-8]	100	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Iron (Fe)	[7439-89-6]	500	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Lead (Pb)	[7439-92-1]	50	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Manganese (Mn)	[7439-96-5]	100	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Nickel (Ni)	[7440-02-0]	100	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Selenium (Se)	[7782-49-2]	25	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Silver (Ag)	[7440-22-4]	25	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Thallium (Tl)	[7440-28-0]	25	µg/mL		
ACCU-ICP-MS-SPIKE-W-1	Vanadium (V)	[7440-62-2]	100	µg/mL		

ACCU-ICP-MS-SPIKE-W-1 Zinc (Zn)

[7440-66-6]

250

µg/mL