

ICP-MS Tuning Solution, 100 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-ICP-MS-TUNSOL1-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	100	µg/mL		
	Beryllium (Be)	[7440-41-7]	100	µg/mL		
	Copper (Cu)	[7440-50-8]	100	µg/mL		
	Indium (In)	[7440-74-6]	100	µg/mL		
	Lithium (Li)	[7439-93-2]	100	µg/mL		
	Magnesium (Mg)	[7439-95-4]	100	µg/mL		
	Thallium (Tl)	[7440-28-0]	100	µg/mL		
	Uranium (U)	[7440-61-1]	100	µg/mL		