

## **Quality Control Standard 23, 1000 µg/mL in 2-5% Nitric Acid**

Art. ID ACCU-JY-QC23-ASL-1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Silver (Ag)	[7440-22-4]	1000	µg/mL		
	Aluminium (Al)	[7429-90-5]	1000	µg/mL		
	Boron (B)	[7440-42-8]	1000	µg/mL		
	Barium (Ba)	[7440-39-3]	1000	µg/mL		
	Bismuth (Bi)	[7440-69-9]	1000	µg/mL		
	Cadmium (Cd)	[7440-43-9]	1000	µg/mL		
	Calcium (Ca)	[7440-70-2]	1000	µg/mL		
	Chromium (Cr)	[7440-47-3]	1000	µg/mL		
	Cobalt (Co)	[7440-48-4]	1000	µg/mL		
	Copper (Cu)	[7440-50-8]	1000	µg/mL		
	Iron (Fe)	[7439-89-6]	1000	µg/mL		
	Gallium (Ga)	[7440-55-3]	1000	µg/mL		
	Indium (In)	[7440-74-6]	1000	µg/mL		
	Potassium (K)	[7440-09-7]	1000	µg/mL		
	Lithium (Li)	[7439-93-2]	1000	µg/mL		
	Magnesium (Mg)	[7439-95-4]	1000	µg/mL		
	Manganese (Mn)	[7439-96-5]	1000	µg/mL		
	Sodium (Na)	[7440-23-5]	1000	µg/mL		
	Nickel (Ni)	[7440-02-0]	1000	µg/mL		
	Lead (Pb)	[7439-92-1]	1000	µg/mL		
	Strontium (Sr)	[7440-24-6]	1000	µg/mL		
	Thallium (Tl)	[7440-28-0]	1000	µg/mL		
	Zinc (Zn)	[7440-66-6]	1000	µg/mL		