

Quality Control Standard 23, 1000 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-JY-QC23-ASL-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Silver (Ag) | [7440-22-4] | 1000 | µg/mL | | |
| | Aluminium (Al) | [7429-90-5] | 1000 | µg/mL | | |
| | Boron (B) | [7440-42-8] | 1000 | µg/mL | | |
| | Barium (Ba) | [7440-39-3] | 1000 | µg/mL | | |
| | Bismuth (Bi) | [7440-69-9] | 1000 | µg/mL | | |
| | Cadmium (Cd) | [7440-43-9] | 1000 | µg/mL | | |
| | Calcium (Ca) | [7440-70-2] | 1000 | µg/mL | | |
| | Chromium (Cr) | [7440-47-3] | 1000 | µg/mL | | |
| | Cobalt (Co) | [7440-48-4] | 1000 | µg/mL | | |
| | Copper (Cu) | [7440-50-8] | 1000 | µg/mL | | |
| | Iron (Fe) | [7439-89-6] | 1000 | µg/mL | | |
| | Gallium (Ga) | [7440-55-3] | 1000 | µg/mL | | |
| | Indium (In) | [7440-74-6] | 1000 | µg/mL | | |
| | Potassium (K) | [7440-09-7] | 1000 | µg/mL | | |
| | Lithium (Li) | [7439-93-2] | 1000 | µg/mL | | |
| | Magnesium (Mg) | [7439-95-4] | 1000 | µg/mL | | |
| | Manganese (Mn) | [7439-96-5] | 1000 | µg/mL | | |
| | Sodium (Na) | [7440-23-5] | 1000 | µg/mL | | |
| | Nickel (Ni) | [7440-02-0] | 1000 | µg/mL | | |
| | Lead (Pb) | [7439-92-1] | 1000 | µg/mL | | |
| | Strontium (Sr) | [7440-24-6] | 1000 | µg/mL | | |
| | Thallium (Tl) | [7440-28-0] | 1000 | µg/mL | | |
| | Zinc (Zn) | [7440-66-6] | 1000 | µg/mL | | |