

## **Density Standard 0.7033g/ml @20°C**

|                 |  |
|-----------------|--|
| Art. ID         | REA-DEN20020   |
| Unit            | 100 mL   |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2 |

### Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 20°C          | Density           |            | 0,7033              | g/mL |        |        |