

KDWR VOC Mix A, 100 µg/mL in Methanol

Art. ID ACCU-KDWR-001-PAK
Unit 5 x 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene	[71-43-2]	100	µg/mL		
	Bromodichloromethane	[75-27-4]	100	µg/mL		
	Bromoform (Tribromomethane)	[75-25-2]	100	µg/mL		
	Chloroform (Trichloromethane)	[67-66-3]	100	µg/mL		
	Dibromochloromethane	[124-48-1]	100	µg/mL		
	Ethylbenzene	[100-41-4]	100	µg/mL		
	Methylene chloride (Dichloromethane)	[75-09-2]	100	µg/mL		
	Phenol	[108-95-2]	100	µg/mL		
	Tetrachloroethene	[127-18-4]	100	µg/mL		
	Toluene	[108-88-3]	100	µg/mL		
	1,1,1-Trichloroethane	[71-55-6]	100	µg/mL		
	Trichloroethene	[79-01-6]	100	µg/mL		
	m-Xylene	[108-38-3]	100	µg/mL		
	p-Xylene	[106-42-3]	100	µg/mL		
	o-Xylene	[95-47-6]	100	µg/mL		