

Density Standard 1.0005g/ml @20°C

| | |
|-----------------|---|
| Art. ID | REA-DEN20140 |
| Unit | 100 mL |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1105 III E1 |

Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 20°C | Density | | 1,0005 | g/mL | | |