

Method 200.7 Instrument Fortifying Solution 2, 20 µg/mL in 2-5% Nitric Acid, tr Hydrofluoric acid

Art. ID ACCU-M-200.7-LFSS-02-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	20	µg/mL		
	Molybdenum (Mo)	[7439-98-7]	20	µg/mL		
	Silicon (Si)	[7440-21-3]	20	µg/mL		
	Tin (Sn)	[7440-31-5]	20	µg/mL		
	Titanium (Ti)	[7440-32-6]	20	µg/mL		