

**Spiking Std #2 Method 200.7 (1991), 10000 µg/mL in 2-5% Nitric Acid**

Art. ID ACCU-M-200.7-SP-02-R  
Unit 50 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Calcium (Ca)	[7440-70-2]	10000	µg/mL		
	Magnesium (Mg)	[7439-95-4]	10000	µg/mL		
	Potassium (K)	[7440-09-7]	10000	µg/mL		
	Sodium (Na)	[7440-23-5]	10000	µg/mL		