

Density Standard 0.6872g/ml @40°C

| | |
|-----------------|------------------------------------------------------------------|
| Art. ID | REA-DEN40020 |
| Unit | 100 mL |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1262 II E2 |

Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 40°C | Density | | 0,6872 | g/mL | | |