

### **EPA Method 524,3 Internal Standard, 2,0 mg/mL in Methanol**

Art. ID ACCU-M-524R-C-IS  
Unit 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,4-Difluorobenzene	[540-36-3]	2	mg/mL		
	Chlorobenzene D5	[3114-55-4]	2	mg/mL		
	1,4-Dichlorobenzene D4	[3855-82-1]	2	mg/mL		