

### **EPA Method 524,3 Surrogate Standard, 2,0 mg/mL in Methanol**

Art. ID ACCU-M-524R-C-SS-PAK  
Unit 5 x 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter                        | CAS number   | Concentration/Value | Unit  | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
|                  | tert.-Butyl-methylether<br>-D3 (MTBE-D3) | [29366-08-3] | 2                   | mg/mL |        |        |
|                  | 1,2-Dichlorobenzene D4                   | [2199-69-1]  | 2                   | mg/mL |        |        |
|                  | 4-Bromofluorobenzene                     | [460-00-4]   | 2                   | mg/mL |        |        |