

EPA Method 524,3 Surrogate Standard, 2,0 mg/mL in Methanol

Art. ID ACCU-M-524R-C-SS-PAK
Unit 5 x 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|--|--------------|---------------------|-------|--------|--------|
| | tert.-Butyl-methylether -D3 (MTBE-D3) | [29366-08-3] | 2 | mg/mL | | |
| | 1,2-Dichlorobenzene D4 | [2199-69-1] | 2 | mg/mL | | |
| | 4-Bromofluorobenzene | [460-00-4] | 2 | mg/mL | | |