

### **EPA Method 524,3 Surrogate Standard, 2,0 mg/mL in Methanol**

Art. ID ACCU-M-524R-C-SS-PAK  
Unit 5 x 1 mL  
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	tert.-Butyl-methylether -D3 (MTBE-D3)	[29366-08-3]	2	mg/mL		
	1,2-Dichlorobenzene D4	[2199-69-1]	2	mg/mL		
	4-Bromofluorobenzene	[460-00-4]	2	mg/mL		