

## **Method 525.2 - Surrogate Standard, 500 µg/mL in Acetone**

Art. ID ACCU-M-525.2-SS2-ASL-PAK  
Unit 5 x 1 mL  
Deliverydetails De Minimis UN1090 II E2

| Text/Information | Analyte/Parameter           | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-----------------------------|-------------|---------------------|-------|--------|--------|
|                  | 1,3-Dimethyl-2-nitrobenzene | [81-20-9]   | 500                 | µg/mL |        |        |
|                  | Perylene D12                | [1520-96-3] | 500                 | µg/mL |        |        |
|                  | Pyrene D10                  | [1718-52-1] | 500                 | µg/mL |        |        |
|                  | Triphenylphosphate          | [115-86-6]  | 500                 | µg/mL |        |        |