

Method 525 - Internal Standard, 2,0 mg/mL in Methylene chloride (Dichloromethane)

Art. ID ACCU-M-525-IS-PAK
Unit 5 x 1 mL
Deliverydetails De Minimis UN1593 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene D10	[15067-26-2]	2	mg/mL		
	Chrysene D12	[1719-03-5]	2	mg/mL		
	Perylene D12	[1520-96-3]	2	mg/mL		
	Phenanthrene D10	[1517-22-2]	2	mg/mL		