

Method 529 Surrogate Analyte Fortification Solution, 100 µg/mL in Methanol

Art. ID ACCU-M-529-SAFS
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,3,5-Trimethyl-2-nitro benzene	[603-71-4]	100	µg/mL		
	1,2,4-Trimethyl-5-nitro benzene	[610-91-3]	100	µg/mL		
	Nitrobenzene-d5	[4165-60-0]	100	µg/mL		