

## **Density Standard 0.7454g/ml @50°C**

|                 |  |
|-----------------|--|
| Art. ID         | REA-DEN50070   |
| Unit            | 100 mL   |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1294 II E2 |

### Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 50°C          | Density           |            | 0,7454              | g/mL |        |        |