

**o,p-DDD & o,p-DDE & o,p-DDT, 0,25 mg/mL in Hexane:Toluene (1:1)**

Art. ID ACCU-M-8080-OP  
Unit 1 mL  
Deliverydetails De Minimis UN1993 II E2

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
|                  | 2,4'-DDD          | [53-19-0]   | 250                 | µg/mL |        |        |
|                  | 2,4'-DDE          | [3424-82-6] | 250                 | µg/mL |        |        |
|                  | 2,4'-DDT          | [789-02-6]  | 250                 | µg/mL |        |        |