

### **Density Standard 0.8098g/ml @80°C**

Art. ID	REA-DEN80110
Unit	100 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1294 II E2

#### Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
at 80°C	Density		0,8098	g/mL		