

## **Density Standard 0.9006g/ml @80°C**

|                 |   |
|-----------------|---|
| Art. ID         | REA-DEN80130  |
| Unit            | 100 mL  |
| Deliverydetails | ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1105 III E1 |

### Description

Free from Toxic heavy Metals in formulation. Manufactured to meet ASTM protocols. Tested, certified and batch traceable. Current Certificates are 17025 approved.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
| at 80°C          | Density           |            | 0,9006              | g/mL |        |        |