

Method 8270 - Surrogate Standard, Varied conc, in Dichloromethane:Methanol (80:20)

Art. ID ACCU-M-8270-SS-R-PAK
Unit 5 x 1 mL
Deliverydetails De Minimis UN1992 (6.1) II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Fluorobiphenyl	[321-60-8]	1000	µg/mL		
	2-Fluorophenol	[367-12-4]	2000	µg/mL		
	Nitrobenzene D5	[4165-60-0]	1000	µg/mL		
	Phenol D5	[4165-62-2]	2000	µg/mL		
	p-Terphenyl D14	[1718-51-0]	1000	µg/mL		
	2,4,6-Tribromophenol	[118-79-6]	2000	µg/mL		