

Multi-Method Calibration Std #1, Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-MCS-01-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Beryllium (Be)	[7440-41-7]	50	µg/mL		
	Cadmium (Cd)	[7440-43-9]	150	µg/mL		
	Lead (Pb)	[7439-92-1]	500	µg/mL		
	Manganese (Mn)	[7439-96-5]	100	µg/mL		
	Selenium (Se)	[7782-49-2]	200	µg/mL		
	Zinc (Zn)	[7440-66-6]	150	µg/mL		