

Multi-Method Calibration Std #2, Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-MCS-02-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	100	µg/mL		
	Cobalt (Co)	[7440-48-4]	100	µg/mL		
	Copper (Cu)	[7440-50-8]	100	µg/mL		
	Iron (Fe)	[7439-89-6]	10000	µg/mL		
	Vanadium (V)	[7440-62-2]	100	µg/mL		