

### **Multi-Method Calibration Std #4R, Varied conc, in 2-5% Nitric Acid**

Art. ID ACCU-MCS-04R-5  
Unit 500 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	200	µg/mL		
	Calcium (Ca)	[7440-70-2]	1000	µg/mL		
	Chromium (Cr)	[7440-47-3]	20	µg/mL		
	Lithium (Li)	[7439-93-2]	100	µg/mL		
	Nickel (Ni)	[7440-02-0]	20	µg/mL		
	Potassium (K)	[7440-09-7]	400	µg/mL		
	Sodium (Na)	[7440-23-5]	200	µg/mL		
	Strontium (Sr)	[7440-24-6]	10	µg/mL		