

Multi-Method Calibration Std #5R, Varied conc, in 2-5% Nitric Acid

Art. ID ACCU-MCS-05R-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Antimony (Sb)	[7440-36-0]	200	µg/mL		
	Magnesium (Mg)	[7439-95-4]	1000	µg/mL		
	Silver (Ag)	[7440-22-4]	50	µg/mL		
	Thallium (Tl)	[7440-28-0]	200	µg/mL		