

Mixed Calibration Standard 6R, Varied conc, in 2-5% Nitric Acid, tr Hydrofluoric acid

Art. ID ACCU-MCS-06R-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Phosphorus (P)	[7723-14-0]	200	µg/mL		
	Tin (Sn)	[7440-31-5]	200	µg/mL		
	Titanium (Ti)	[7440-32-6]	100	µg/mL		
	Boron (B)	[7440-42-8]	50	µg/mL		
	Silicon (Si)	[7440-21-3]	100	µg/mL		