

ICP multi-element standard solution IX, 100 µg/mL in 1 mol/L Nitric Acid

Art. ID ACCU-MES-09-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN2031 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Arsenic (As)	[7440-38-2]	100	µg/mL		
	Beryllium (Be)	[7440-41-7]	100	µg/mL		
	Lead (Pb)	[7439-92-1]	100	µg/mL		
	Cadmium (Cd)	[7440-43-9]	100	µg/mL		
	Chromium (Cr)	[7440-47-3]	100	µg/mL		
	Nickel (Ni)	[7440-02-0]	100	µg/mL		
	Selenium (Se)	[7782-49-2]	100	µg/mL		
	Thallium (Tl)	[7440-28-0]	100	µg/mL		