

ASTM/EPA Gasoline Refinery Internal Standard Mix, Varied Conc. by Weight

Art. ID ACCU-M-GRA-IS-5ML-PAK
Unit 5 x 5 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) to be checked

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene D6	[1076-43-3]	40	%		
	Ethylbenzene-D10	[25837-05-2]	40	%		
	Naphthalene D8	[1146-65-2]	20	%		