

ASTM/EPA Gasoline Refinery Internal Standard Mix, Varied conc, by %w/w

Art. ID ACCU-M-GRA-IS-R-10ML
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) to be checked

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Benzene D6	[1076-43-3]	16,67	g/100g		
	Ethylbenzene-D10	[25837-05-2]	16,65	g/100g		
	Naphthalene D8	[1146-65-2]	8,77	g/100g		
	Toluene D8	[2037-26-5]	57,91	g/100g		