

Instrument Calibration Standard 3, 1000 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-PE-CAL3-ASL-5
Unit 500 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Iron (Fe)	[7439-89-6]	1000	µg/mL		
	Potassium (K)	[7440-09-7]	1000	µg/mL		
	Calcium (Ca)	[7440-70-2]	1000	µg/mL		
	Sodium (Na)	[7440-23-5]	1000	µg/mL		
	Magnesium (Mg)	[7439-95-4]	1000	µg/mL		