

Multi-element Calibration Standard 1, 10 µg/mL in 2-5% Nitric Acid

Art. ID ACCU-PE-MECAL1-ASL-1
Unit 100 mL
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
| | Beryllium (Be) | [7440-41-7] | 10 | µg/mL | | |
| | Bismuth (Bi) | [7440-69-9] | 10 | µg/mL | | |
| | Cerium (Ce) | [7440-45-1] | 10 | µg/mL | | |
| | Cobalt (Co) | [7440-48-4] | 10 | µg/mL | | |
| | Indium (In) | [7440-74-6] | 10 | µg/mL | | |
| | Magnesium (Mg) | [7439-95-4] | 10 | µg/mL | | |
| | Nickel (Ni) | [7440-02-0] | 10 | µg/mL | | |
| | Lead (Pb) | [7439-92-1] | 10 | µg/mL | | |
| | Uranium (U) | [7440-61-1] | 10 | µg/mL | | |