

*****No longer available - see new batch NW-TMDA-64.3***Water - Trace elements**

Art. ID NW-TMDA-64.2
Unit 250 mL
Deliverydetails No Dangerous Good /not restricted

Description
lot 0313

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	290 ± 26,5	µg/L		
	Antimony (Sb)	[7440-36-0]	128 ± 12,2	µg/L		
	Arsenic (As)	[7440-38-2]	162 ± 15,4	µg/L		
	Barium (Ba)	[7440-39-3]	290 ± 23,5	µg/L		
	Beryllium (Be)	[7440-41-7]	160 ± 15,6	µg/L		
	Bismuth (Bi)	[7440-69-9]	133 ± 27,3	µg/L		
	Boron (B)	[7440-42-8]	282 ± 23,8	µg/L		
	Cadmium (Cd)	[7440-43-9]	266 ± 21,3	µg/L		
	Chromium (Cr)	[7440-47-3]	290 ± 19,6	µg/L		
	Cobalt (Co)	[7440-48-4]	254 ± 17,6	µg/L		
	Copper (Cu)	[7440-50-8]	275 ± 23,3	µg/L		
For Information	Gallium (Ga)	[7440-55-3]	53*	µg/L		
	Iron (Fe)	[7439-89-6]	306 ± 27,4	µg/L		
	Lead (Pb)	[7439-92-1]	288 ± 25,2	µg/L		
	Lithium (Li)	[7439-93-2]	151 ± 16,2	µg/L		
	Manganese (Mn)	[7439-96-5]	295 ± 21,5	µg/L		
	Molybdenum (Mo)	[7439-98-7]	290 ± 25,4	µg/L		
	Nickel (Ni)	[7440-02-0]	263 ± 19,3	µg/L		
For Information	Rubidium (Rb)	[7440-17-7]	30,6 ± 2,9	µg/L		
	Selenium (Se)	[7782-49-2]	155 ± 18,3	µg/L		
For Information	Silver (Ag)	[7440-22-4]	13*	µg/L		
	Strontium (Sr)	[7440-24-6]	641 ± 46,5	µg/L		
	Thallium (Tl)	[7440-28-0]	146 ± 13,8	µg/L		
	Tin (Sn)	[7440-31-5]	289 ± 37	µg/L		
	Titanium (Ti)	[7440-32-6]	128 ± 9,1	µg/L		
For Information	Tungsten (W)	[7440-33-7]	0,09*	µg/L		
	Uranium (U)	[7440-61-1]	142 ± 13,8	µg/L		
	Vanadium (V)	[7440-62-2]	289 ± 21	µg/L		
	Zinc (Zn)	[7440-66-6]	310 ± 27,7	µg/L		

