

SmartTune Solution for DRC / DRCPlus / DRC II, Varied conc, in 0,5% Nitric Acid

Art. ID ACCU-PE-SMTUNE2-ASL-5
Unit 500mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	10	µg/mL		
	Beryllium (Be)	[7440-41-7]	1	µg/mL		
	Cerium (Ce)	[7440-45-1]	1	µg/mL		
	Cobalt (Co)	[7440-48-4]	1	µg/mL		
	Indium (In)	[7440-74-6]	1	µg/mL		
	Iron (Fe)	[7439-89-6]	1	µg/mL		
	Lead (Pb)	[7439-92-1]	1	µg/mL		
	Magnesium (Mg)	[7439-95-4]	1	µg/mL		
	Thorium (Th)	[7440-29-1]	1	µg/mL		
	Uranium (U)	[7440-61-1]	1	µg/mL		