

## **Tuning Solution I, 10 µg/mL in 2-5% Nitric Acid, tr Hydrochloric acid**

Art. ID ACCU-PE-TUNSOL-ASL-1  
Unit 100 mL  
Deliverydetails ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN3264 III E1

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Barium (Ba)	[7440-39-3]	10	µg/mL		
	Beryllium (Be)	[7440-41-7]	10	µg/mL		
	Cerium (Ce)	[7440-45-1]	10	µg/mL		
	Cobalt (Co)	[7440-48-4]	10	µg/mL		
	Indium (In)	[7440-74-6]	10	µg/mL		
	Lithium (Li)	[7439-93-2]	10	µg/mL		
	Magnesium (Mg)	[7439-95-4]	10	µg/mL		
	Lead (Pb)	[7439-92-1]	10	µg/mL		
	Rhodium (Rh)	[7440-16-6]	10	µg/mL		
	Thallium (Tl)	[7440-28-0]	10	µg/mL		
	Uranium (U)	[7440-61-1]	10	µg/mL		
	Yttrium (Y)	[7440-65-5]	10	µg/mL		