

Diethyl Phthalate-3,4,5,6-d4 Standard, 100 µg/mL in Methanol

Art. ID ACCU-PHTH-D4-005S
Unit 1 mL
Deliverydetails De Minimis UN1230 (6.1) II E2

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------------------|--|--------------|---------------------|-------|--------|--------|
| Diethyl Phthalate-3,4,5,6-d4 | Phthalic acid, bis-ethyl ester D4 (Bis-ethyl phthalate-D4, Diethyl phthalate-D4) | [93952-12-6] | 100 | µg/mL | | |