

(S)-(+)-2-Amino-4-bromobutyric Acid Hydrobromide (~90%)

Art. ID TRC-A602465-10G
Unit 10 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Application Note: (S)-(+)-2-Amino-4-bromobutyric Acid Hydrobromide is a building block used in stereoselective preparation of phosphinic acid phosphapeptide analogs of glutamyl-gamma-glutamate as potent inhibitors of folypolyglutamate synthetase. Also used in the synthesis of isoindoline carbamate and acylsulfonamide containing macrocyclic peptidomimetics as inhibitors of HCV NS3/4A protease.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	(S)-(+)-2-Amino-4-bromo butyric acid hydrobromi de	[15159-65-6]				