

4-[(Aminocarbonyl)amino]-D-phenylalanine

Art. ID TRC-A603178-25MG
Unit 25 mg
Deliverydetails De Minimis

Description

Application Note: 4-[(Aminocarbonyl)amino]-D-phenylalanine can be used a reactant in continuous flow solid-phase peptide synthesis.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	4-[(Aminocarbonyl)amino]-D-phenylalanine	[1808153-79-8]				